



CGM Comparison Chart

| Feature | Dexcom G7 ¹ | FreeStyle Libre 3 Plus ² | FreeStyle Libre 2 Plus ³ |
|---|---|-------------------------------------|---|
| Sensor Duration | 10 days and 15 days* | 15 days | 15 days |
| MARD | 8.2% | 8.2% | 8.2% |
| Transmitter Range | 33' (Bluetooth) | 33' (Bluetooth) | 20' (Bluetooth) |
| Warm-Up Time | < 30 min | 1 hour | 1 hour |
| Insulin Pump Integration | <ul style="list-style-type: none"> • iLet Bionic Pancreas* • Lilly Tempo Insulin Pen • Omnipod® 5* • Tandem Mobi System • Tandem t-slim:X2™ insulin pump | iLet Bionic Pancreas | <ul style="list-style-type: none"> • Omnipod® 5 Automated Insulin Delivery System • Tandem t-slim:X2™ |
| Smart Phone Compatible[†] | iPhone and Android | iPhone and Android | iPhone and Android |
| Water Resistance | 2.4 meters for 24 hours | 3 feet for 30 minutes | 3 feet for 30 minutes |
| Approved Age | 2 years and older | 2 years and older | 2 years and older |
| Insulin Pump Integration | Yes | Yes | Yes |

¹ Dexcom. (2025, March 3). Dexcom G7 vs. FreeStyle Libre 3 and Libre 3 Plus: Choosing the Right CGM. Retrieved from <https://www.dexcom.com/all-access/dexcom-cgm-explained/dexcom-g7-vs-freestyle-libre-3-and-libre-3-plus-choosing-the-right-cgm>

² Abbott. (n.d.). FreeStyle Libre 3 Continuous Glucose Monitoring. Retrieved April 28, 2025, from <https://www.freestyleprovider.abbott/us-en/freestyle-libre-3.html>

³ Abbott. (n.d.). FreeStyle Libre 2 Continuous Glucose Monitoring. Retrieved April 28, 2025, from <https://www.freestyleprovider.abbott/us-en/freestyle-libre-2.html>

* For complete product features, please visit the manufacturers website.

[†] Smartphone compatibility may vary by device.

* Dexcom G7 15-day is currently only available for integration with iLet Bionic Pancreas and Omnipod 5.

Better Service, Better Care for Your Patients. Become a Prescriber Partner Today.

Visit www.hcp.advanceddiabetes.com