



CGM Comparison Chart

Feature	Dexcom G7 ¹	FreeStyle Libre 3 Plus ²	FreeStyle Libre 2 Plus ³
Sensor Duration	10 days and 15 days [*]	15 days	15 days
MARD	8.2%	8.2%	8.2%
Transmitter Range	33' (Bluetooth)	33' (Bluetooth)	20' (Bluetooth)
Warm-Up Time	< 30 min	1 hour	1 hour
Insulin Pump Integration	<ul style="list-style-type: none"> • iLet Bionic Pancreas[*] • Lilly Tempo Insulin Pen • Omnipod[®] 5[*] • Tandem Mobi System • Tandem t-slim:X2[™] insulin pump 	iLet Bionic Pancreas	<ul style="list-style-type: none"> • Omnipod[®] 5 Automated Insulin Delivery System • Tandem t-slim:X2[™]
Smart Phone Compatible[†]	iPhone and Android	iPhone and Android	iPhone and Android
Water Resistance	2.4 meters for 24 hours	3 feet for 30 minutes	3 feet for 30 minutes
Approved Age	2 years and older	2 years and older	2 years and older
Insulin Pump Integration	Yes	Yes	Yes

¹ Dexcom. (2025, March 3). Dexcom G7 vs. FreeStyle Libre 3 and Libre 3 Plus: Choosing the Right CGM. Retrieved from <https://www.dexcom.com/all-access/dexcom-cgm-explained/dexcom-g7-vs-freestyle-libre-3-and-libre-3-plus-choosing-the-right-cgm>

² Abbott. (n.d.). FreeStyle Libre 3 Continuous Glucose Monitoring. Retrieved April 28, 2025, from <https://www.freestyleprovider.abbott/us-en/freestyle-libre-3.html>

³ Abbott. (n.d.). FreeStyle Libre 2 Continuous Glucose Monitoring. Retrieved April 28, 2025, from <https://www.freestyleprovider.abbott/us-en/freestyle-libre-2.html>

^{*} For complete product features, please visit the manufacturers website.

[†] Smartphone compatibility may vary by device.

^{*} Dexcom G7 15-day is currently only available for integration with Tandem t-slim:X2, iLet Bionic Pancreas and Omnipod 5.

Better Service, Better Care for Your Patients. Become a Prescriber Partner Today.

Visit hcp.advanceddiabetessupply.com