





Insulin Pump Comparison Chart

Choosing the right insulin pump for your patient may be challenging. To simplify the process, we've created a guide that provides essential information, that may help you find the best fit to support your doctor-patient relationship.




Insulet Omnipod® 5



Insulet Omnipod® Dash



MiniMed™ 780G Pump System



Tandem T-slim X2



Tandem Mobi



iLet Bionic Pancreas

Features and Compatibility

This table outlines key features and compatibility details that may help prescribers quickly compare insulin pump options and choose the best fit for each patient.

Pump Name	Delivery System	Infusion Set (If Tubed)	Reservoir/ Cartridge	Water Resistance
Insulet Omnipod® Dash*¹	Tubeless	Not Needed	200 unit Pod	Pods are waterproof up to 25 ft for 60 min; PDM controller is not waterproof
Insulet Omnipod® 5*¹	Tubeless	Not Needed	200 unit Pod	Pods are waterproof up to 25 ft for 60 min; PDM controller is not waterproof
Tandem T- slim X2²	Tubed	Tandem TruSteel VariSoft, Auto Soft (30, 90, XC) 5" tubing not compatible with X2.	300 unit Cartridge	Watertight Up to 3 ft for 30 min
Tandem Mobi²	Tubed (on-body)	Tandem AutoSoft 30, 60, XC (5" tubing)	200 unit Cartridge	Up to 8 ft for 2 hours
MiniMed™ 780G³	Tubed	Medtronic Quick-set, Sure-T Silhouette, Mio Advance Medtronic Extended Infusion Set (EIS)	<ul style="list-style-type: none"> • 180 unit Reservoir • 300 unit Reservoir • 300 unit Extended Reservoir (EIS) 	Water proof up to 12 ft for 24 hours
iLet Bionic Pancreas⁴	Tubed	Inset, Contact Detach	180 unit Cartridge	Protected against immersion in water up to 12 ft for 30 min

*Because of the nature of the Omnipod systems (sold as Kit + pods), you can simply submit your patient's prescription. There are no additional accessories required.

In case of pump or CGM failure, patients will need back up long-acting insulin, syringes/pen needles, and blood glucose meter.

Helpful Tip: Choose infusion sets which align with patient's body type, activity level and personal comfort preferences, or add "patient preference" for model and size on the prescription.

Better Service, Better Care for Your Patients. Become a Prescriber Partner Today.

Visit hcp.advanceddiabetessupply.com

References:

1. Insulet Corporation (www.omnipod.com)
2. Tandem Diabetes (www.tandemdiabetes.com)
3. Medtronic Diabetes (www.medtronicdiabetes.com)
4. Beta Bionics (www.betabionics.com)